

Report Number: 030092
Revision Number: 1

Page 1 of 7

**Determination of the Antimicrobial Activity
Of A Disinfectant Solution Based on BS 1276:1997**

Test Product: UB20 Sanitiser Concentrate
Silliker Reference No. 03040171

Client Name: Des Eustace
Clover Chemicals
Clover House
Macclesfield Road
Whaley Bridge
High Peak
SK23 7DQ

Date Received: 7th April 2003

Date Testing Commenced: 11th April 2003

Date Completed: 25th April 2003

Author: D Prince

Authorised By:.....Date:.....
Mr Carl Winfield BSc
Microbiology Laboratory Manager

Circulation List

Original: Clover Chemicals
Copy Number 1: Silliker Archives

SILLIKER LTD.
PARK ROAD
OVERSEAL
SWADLINCOTE
DE12 6JT
TEL. 01283 213036
FAX. 01283 550713

REGISTERED IN ENGLAND No. 3624497

1. **Scope**

The following report gives details of the method used to determine the antimicrobial activity of disinfectants used in Food Hygiene and is based on British Standard EN 1276:1997.

2. **Test Products**

2.1 UB20 Sanitiser Concentrate Silliker Reference No. 03040171

3. **Test Organisms**

3.1	<i>Staphylococcus aureus</i>	NCTC 10788	Silliker Ref S3
3.2	<i>Enterococcus hirae</i>	ATCC 8043	Silliker Ref E10
3.3	<i>Pseudomonas aeruginosa</i>	NCIMB 10421	Silliker Ref P2
3.4	<i>Escherichia coli</i>	ATCC 10536	Silliker Ref E9

4. **Equipment / Materials**

- 4.1 Diluent
- 4.2 Bovine albumin solution 0.3g/100ml
- 4.3 Water of Standard Hardness (WSH)
- 4.4 Inactivation Liquid
- 4.5 Tryptone Soya Agar (TSA)

All above prepared as per EN 1276:1997

- 4.6 Incubator at 36°C +/- 1°C
- 4.7 Sterile 1ml Pipettes
- 4.8 Sterile 10ml Pipettes
- 4.9 Sterile 90mm Petri Dishes
- 4.10 Sterile 30ml plastic
- 4.11 Incubator at 37°C ± 1°C

5. **Preparation of Test Microbial Suspensions**

- 5.1 For each organism, inoculate the surface of Tryptone Soya Agar plates from a stock culture. Incubate the plates at 37°C +/- 1°C for 18 – 24 hours.
- 5.2 Make a suspension containing approximately 1.5×10^8 to 5×10^8 cfu/ml for each micro-organism. Serially dilute the suspension down in diluent to obtain a count of between 6.0×10^2 and 3×10^3 cfu/ml. Pour with TSA and incubate at 36°C +/- 1°C for 48 hours. This is used as a control count.

6. Preparation of Test Disinfectant Solutions

A volumetric solution 1.25 times stronger than the required concentration of test product was prepared in WSH.

7. Determination of the Microbial Effect of the Disinfectant Solutions

7.1 Clean Conditions (Organic Load 0.03% Bovine Albumin)

7.1.1 1ml of 0.3g/100ml bovine albumin solution was added to each clearly labelled sterile 30ml universal. 1ml of the appropriate microbial suspension was added to each container and mixed. After 2 minutes incubation at $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$. 8mls of product solution was added to each container and mixed. After 5 minutes ± 10 seconds (contact time) at $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$, 1ml of liquid from each container was transferred into separate 8ml aliquots of inactivation liquid and 1ml of WSH and mixed. After a neutralisation time of 5 minutes ± 10 s, 1ml of the test solution was serially diluted out in diluent and plated out into duplicate petri dishes.

Approximately 15mls of TSA melted and tempered to $45^{\circ}\text{C} \pm 1^{\circ}\text{C}$ was added to the plates, mixed and allowed to set before incubating at $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 48 hours.

7.1.2 This procedure was repeated for each of the test microbial suspensions. And for each product dilution.

8. Assay Validation

To ensure that the experimental conditions do not interfere with the test, the following was done:

Temperature Validation

A suspension containing between 6×10^2 and 3×10^3 of each test organism was prepared in diluent. 1ml of the interfering substance and 1ml of the bacterial suspension was added to a sterile universal and mixed. After an incubation time of 2 min ± 10 s at $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$. 8.0mls of hard water was added and the suspension mixed. After an incubation time of 5 min ± 10 s at $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$, the test solution was plated out into duplicate petri dishes. The above procedure was carried out for each organism at the highest product concentration. The plates were poured with 15ml of TSA and incubated at $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$.

Neutraliser Validation

A suspension containing between 6×10^2 and 3×10^3 of each test organism was prepared in diluent. 1ml of the bacterial suspension was added to a universal already containing 8ml of neutraliser and 1ml of hard water and mixed. After a contact time of 5 min ± 10 s at $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$, 1ml of the test solution was plated out into duplicate petri dishes. The above procedure was carried out for each organism at the highest product concentration. The plates were poured with 15ml of TSA and incubated at $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$.

Dilution / Neutralisation Validation

1ml of the interfering substance and 1ml of diluent was added to a sterile universal and mixed. After starting a stopwatch, 1ml of the product dilution was added and mixed. After a contact time of 5 min \pm 10s at 20°C \pm 1°C, 1ml of the test mixture was added to 8ml of neutraliser for 5 min \pm 10s. 1ml of the bacterial suspension containing between 6 x 10² and 3 x 10³ cfu/ml was added to the test mixture and left for a contact time of 30min \pm 10s at 20°C \pm 1°C. After the contact time, 1ml of the test solution was plated out into duplicate petri dishes. The above procedure was carried out for each organism at the highest product concentration. The plates were poured with 15ml of TSA and incubated at 36°C \pm 1°C.

Control Counts

Organism	Dilution	Count on Plate		
		1	2	Mean
P2	-6	112	120	1.2 x 10 ⁸
	-7	14	24	1.9 x 10 ⁸
S3	-6	212	220	2.2 x 10 ⁸
	-7	40	30	3.5 x 10 ⁸
E9	-6	240	276	2.6 x 10 ⁸
	-7	25	36	3.1 x 10 ⁸
E10	-6	220	200	2.1 x 10 ⁸
	-7	30	18	2.4 x 10 ⁸

Validation of Neutraliser

Organism	Count on Plate		
	1	2	Mean
P2	35	29	3.2 x 10 ²
S3	36	32	3.4 x 10 ²
E9	30	30	3.0 x 10 ²
E10	23	31	2.7 x 10 ²

Temperature Validation

Organism	Count on Plate		
	1	2	3
P2	37	41	3.9 x 10 ²
S3	46	35	4.1 x 10 ²
E9	38	45	4.2 x 10 ²
E10	36	30	3.3 x 10 ²

Validation – Dilution / Neutralisation

Sample	Organism	Count on Plate		
		1	2	Mean
0304 0171	P2	10	12	11
0304 0172		9	10	10
0304 0717	S3	8	12	10
0304 0172		9	8	9
0304 0171	E9	25	20	23
0304 0172		20	22	21
0304 0171	E10	15	10	13
0304 0172		22	18	20

UB20 Sanitiser Concentrate Silliker Reference No. 03040171

Organisms	Initial Inoculum	1:75			1:100			1:125			1:150			1:175			1:200		
		1	2	Mean	1	2	Mean	1	2	Mean	1	2	Mean	1	2	Mean	1	2	Mean
P2	1.2 x 10 ⁸	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10
S3	2.2 x 10 ⁸	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10
E9	3.1 x 10 ⁸	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10
E10	2.1 x 10 ⁸	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10	< 10

Report Number: 030092
Revision Number: 1

Page 7 of 7

Conclusion

A 5 log reduction from the initial count must be achieved for the product to pass the set criteria.

UB20 Sanitiser Concentrate Siliker Reference 0304 0171 passed the criteria for all test dilutions and micro-organisms tested.